



MTA micro placement instrument



Testimonials

When you are able to place MTA exactly where it is needed, you finally understand how essential the Map System is.

Dr. Riccardo Tonini

The flexibility of NiTi needles is a real added value. The various MAP System kits are dedicated not only to endodontic specialists, but will also be very useful for general practitioners whomay need to use these devices in their daily practice.

Dr. Stéphane Simon

Thanks to the MAP System, I have the possibility to achieve a perfect retrograde filling of 4 mm deep without overflow.

The quality of the filling is significantly improved. The triple angle needles provide excellent access and a perfect overview of the surgical field.

Dr. Bernd Ilgenstein

I have been using MAP System formany years for MTA placement. Its clinical versatility is incomparable, especially due to the NiTi bendable tips that allow the best ergonomoy in delivering the material in themost challenging scenarios.

Dr. Marco Martignoni

MAP System is simply amazing. It allows clinicians to pre-bend the NiTi carrier and place MTA without any problems since the plugger inside the tube is plastic and smoothly follows the curvature. I find MAP systemremarkable, not only for the intra-canal use, but also for pulp capping procedures that demand precision and confidence. It's simply genius!

Dr. Alexandru Gliga

MAP System

SMART DESIGN FOR MODERN DENTISTRY

The must-have solution for simple and precise MTA placement.

Efficiency

Unique method for precisely placing MTA by orthograde or retrograde approaches. This professional system is smart and considered the product of choice for perforation, root-end filling, pulp-capping and retro-obturation therapies.

Precision

Wide range and sizes of color coded needles to cover all indications. NiTi Memory Shape needles are by design universal and ideal for both orthograde and retrograde therapies while Classic curved needles are designed for very precise non-surgical endodontic procedures. The Surgery triple angle and Hooked needles are best indicated for surgical endodontic procedures.

Comfort

Plastic Plungers are single use and color coded for easy identification. The Paste Cleaning O-rings are color coded for easy identification and delivered with a customized key for an ergonomic and easy removal of excess material around the needle. Finally, the Cleaning Curettes are sterilizable and made of NiTi with a sharp tip to help remove residues of material inside the needles.

MAP System is made in Switzerland, produced by experienced and highly qualified technicians.


- Extremely high-quality level of material and manufacturing process
- Highly durable and cost-effective
- Color coded for simplified diameter identification

NiTi Memory Shape

This memory shape needle can be manually shaped to any required curvature to provide perfect access during orthograde and retrograde procedures. The needle takes back its initial straight shape during sterilization.

Indications

- Direct pulp capping
- Root canal obturation in intersession or during apexogenesis and revascularization procedures
- Apexification of immature teeth
- Repair of iatrogenic or pathological root canal perforations
- Obturation via surgical access of cervical radicular resorptions
- Simple retrograde obturation after apical resection





Surgery

This triple-angle needle enables an easy access and perfect vision of the entire surgical field during complex retrograde obturation cases.

Indications

- Retrograde obturation after apical resection
- Obturation via surgical access of cervical radicular resorptions




Classic

This curved needle enables an extremely precise gesture and placement of the material via orthograde access.

Indications

- Direct pulp capping
- Root canal obturation in intersession or during apexogenesis and revascularization procedures
- Apexification of immature teeth
- Repair of iatrogenic or pathological root canal perforations



Hooked

This hooked needle enables an easy access and perfect vision of the entire surgical field during retrograde obturation cases.

Indications

- Retrograde obturation after apical resection


MAP One

REF 20296


- 1 NiTi Memory Shape Needle N° 1,
- 8 Plastic Plungers N° 1,
- 1 Cleaning Curette N° 1,
- 1 Stainless Steel Syringe




MAP Needles




NiTi Memory Shape			
N° 0	Yellow	0.90 mm	20210
N° 1	Red	1.10 mm	20211
N° 2	Blue	1.30 mm	20212



Classic			
N° 0	Yellow	0.90 mm	20220
N° 1	Red	1.10 mm	20221



Surgery			
N° 0 droite	Yellow	0.90 mm	20225
N° 0 gauche	Yellow	0.90 mm	20226
N° 1 droite	Red	1.10 mm	20227
N° 1 gauche	Red	1.10 mm	20228



Hooked			
N° 0	Yellow	0.90 mm	20263
N° 1	Red	1.10 mm	20264
N° 2	Blue	1.30 mm	20265

MAP Kits

Intro Kit (3-needles)

REF 20285

Optimal MAP System kit for all professionals looking for flexibility and precision

1 NiTi Memory Shape Needle N° 0
1 NiTi Memory Shape Needle N° 1
1 NiTi Memory Shape Needle N° 2
8 Plastic Plungers N° 0
8 Plastic Plungers N° 1
8 Plastic Plungers N° 2
1 Cleaning Curette N° 0
1 Cleaning Curette N° 1
1 Cleaning Curette N° 2
1 Stainless Steel Syringe
1 Paste Cleaning O-rings N° 0
1 Paste Cleaning O-rings N° 1
1 Paste Cleaning O-rings N° 2
1 Stainless Steel Well
1 Aluminium Case 185 × 95 × 30 mm



Universal Kit (5-needles)

REF 20286

Complete kit for frequent users of the MAP System who need maximum access and placement options

1 NiTi Memory Shape Needle N° 0
1 NiTi Memory Shape Needle N° 1
1 NiTi Memory Shape Needle N° 2
1 Classic Needle N° 0
1 Classic Needle N° 1
16 Plastic Plungers N° 0
16 Plastic Plungers N° 1
8 Plastic Plungers N° 2
1 Cleaning Curette N° 0
1 Cleaning Curette N° 1
1 Cleaning Curette N° 2
1 Stainless Steel Syringe
1 Paste Cleaning O-rings N° 0
1 Paste Cleaning O-rings N° 1
1 Paste Cleaning O-rings N° 2
1 Stainless Steel Well
1 Aluminium Case 185 × 145 × 30 mm



MAP Accessories

Plastic plunger

N° 0	Yellow	20200
N° 1	Red	20201
N° 2	Blue	20202

Pack of 16 plunger

Paste cleaning O-rings

N° 0	Yellow	20252
N° 1	Red	20253
N° 2	Blue	20254

Box of 50 O-rings

Cleaning curettes

N° 0	Yellow	20240
N° 1	Red	20241
N° 2	Blue	20243